


## R environment setting


If run successfully, you will launch the following GUI:



When running STABox for the first time, please specify the installation path of the R environment:


R\_HOME: D:\\software\\miniconda3\\envs\\STABox\\lib\\R

R\_USER: D:\\software\\miniconda3\\envs\\STABox\\Lib\\site-packages\\rpy2




# Running STAGATE


Firstly, we need to specify the pre-processed h5ad file path, then select STAGATE method and click “Confirm”.




Next, we need to set the training parameters of the STAGATE model, and then click “Run” to start the training process.




During the model training process, we can monitor the training progress by checking the command line window. The results are stored at `D:\software\STABox\src\stabox\result\STAGATE_output.h5ad`, and we can find this path from the print history from the command window.




Next, we import the h5ad file with running results, and run the three clustering methods: mclust, Louvain, and kmeans. Then, click “Confirm” to start the clustering process.




When running finished, click “Show” to display the results of dimensionality reduction, clustering, and trajectory inference as follows:




We can also display the tissue structure of the slice in 3D space by setting a port number.





## Running STAligner

Firstly, we need to switch to multi-slice running mode, then set the directory path storing h5ad files, and select the h5ad files to be integrated.



Next, we need to set the training parameters of the STAlignermodel, and then click “Run” to start the training process.

